


A NEW DIMENSION
 IN LASER BASED HISTOLOGY
 EXTRA LARGE & ADJUSTABLE



Size

Slide Size

Sample Size

Extra Large

76 x 102 mm (3 inch x 4 inch)

up to 66 x 66 mm

Double Standard

76 x 52 mm (3 inch x 2 inch)


up to 42 x 42 mm

Standard


76 x 26 mm (3 inch x 1 inch)

up to 32 x 20 mm

LARGER SLIDES, LARGER SAMPLES, MORE OPTIONS



- Various sample holder inserts
- Sensors for automated detection of sample holder inserts
- Extra large slides available at LLS
- Customized solutions on request



Benefit from the new dimension in laser based histology,
Upgrade your running system!